

## Screw BIFIXE®

### Compression screw



Indication:

## True cervical femur fracture



### The BIFIXE screw

The *BIFIXE* system allows compression of both fragments for true cervical fractures that are not or little displaced from the neck of the femur. This compression ensures a strengthening of the fractured area. This system does not puncture the external tissues while achieving a 20 mm recoil.

The installation of a *BIFIXE system* is performed only in the operating room using a specific ancillary device.



### TECHNICAL FEATURES

#### Screw neck




- Anodized titanium
- Barrel length: 40 mm
- Barrel diameter: 9.50 mm
- Net diameter: 11.90 mm

#### Screw body



- Anodized titanium
- Compression head screw, threaded over 20 mm and fitted with a pin suitable for sockets
- Screw body length: from 55 mm to 100 mm
- Body diameter: 6.50 mm
- Net diameter: 10.10 mm

#### BIFIXE system



- Anodized titanium
- Screw bodies delivered with screw head
- Compression 20 mm MAX
- System length: 75 mm to 120 mm

(TABLES OF REFERENCES PAGE 2)


Class IIb Medical Device

Before any use, read carefully the instructions on the instructions for use and on the label

FT/TR/02/A – Update : 09/01/2021

## Screw BIFIXE®

### Compression screw



Indication:

**True cervical femur fracture**

#### REFERENCES

##### BIFIX screw neck

References	Length
TR03/1A/2/40	40 mm

Sterilization : Gamma radiation 25 kGy

##### BIFIXE screw body

References	Length
TR03/1A/1/55	55 mm
TR03/1A/1/60	60 mm
TR03/1A/1/65	65 mm
TR03/1A/1/70	70 mm
TR03/1A/1/75	75 mm
TR03/1A/1/80	80 mm
TR03/1A/1/85	85 mm
TR03/1A/1/90	90 mm
TR03/1A/1/95	95 mm
TR03/1A/1/100	100 mm

Sterilization : Gamma radiation 25 kGy

##### BIFIXE System (Assembled Neck and Screw Bodies)

References	Assembled length	Maximum compression
TR03/1A/40/55/75	75mm	20 mm
TR03/1A/40/60/80	80 mm	
TR03/1A/40/65/85	85 mm	
TR03/1A/40/70/90	90 mm	
TR03/1A/40/75/95	95 mm	
TR03/1A/40/80/100	100 mm	
TR03/1A/40/85/105	105 mm	
TR03/1A/40/90/110	110 mm	
TR03/1A/40/95/115	115 mm	
TR03/1A/40/100/120	120 mm	

Sterilization: Gamma radiation 25 kGy

Class IIb Medical Device

Before any use, read carefully the instructions on the instructions for use and on the label

FT/TR/02/A – Update : 09/01/2021